Tropical Journal of Pharmaceutical Research November 2023; 22 (11): 2399-2406 ISSN: 1596-5996 (print); 1596-9827 (electronic) © Pharmacotherapy Group, Faculty of Pharmacy, University of Benin, Benin City, 300001 Nigeria.

> Available online at http://www.tjpr.org http://dx.doi.org/10.4314/tjpr.v22i11.21

Review Article

Disintegration testing for orally disintegrating tablets (ODTs): An overview

Rawan SM Bafail*

Department of Pharmaceutics and Pharmaceutical Technology, College of Pharmacy, Taibah University, Madinah Munawara, Saudi Arabia

*For correspondence: Email: rawanpharma56@yahoo.com

Sent for review: 11 August 2023

Revised accepted: 2 November 2023

Abstract

Prescribing orally disintegrating tablet (ODT) dosage form is a common practice. For several populations, including geriatrics and pediatrics, ODTs are one of the most desired dosage forms. They are solid dosage forms designed to disintegrate or dissolve very quickly in contact with saliva as soon as they are placed in oral cavity of a patient. According to United States Pharmacopeia (USP), every dosage form has to pass some quality control tests to ensure uniformity and performance of this dosage form. For ODTs, disintegration test is fundamental to ascertain time required for tablets to break down and release their contents for dissolution and absorption. Additionally, it is a crucial predictive test for determining in vitro/in vivo association. There are no mandatory universal disintegration testing standards for ODTs, however, as stated in the USP. Instead, recent USP refers to individual manufacturer monograph specifications which may vary from one monograph to another. This article discusses several developed disintegration tests and methods for ODTs such as OD-mate testing, basket rack assembly, CCD camera, texture analyzer (TA), special disintegration equipment, prototype disintegration tester (PDT), simple approach, and modified wetting test and elaborates on their advantages.

Keywords: Orally disintegrating tablets (ODTs), USP physical test, Disintegration test, super disintegrants

This is an Open Access article that uses a funding model which does not charge readers or their institutions for access and distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/4.0) and the Budapest Open Access Initiative (http://www.budapestopenaccessinitiative.org/read), which permit unrestricted use, distribution, and reproduction in any medium, provided the original work is properly credited.

Tropical Journal of Pharmaceutical Research is indexed by Science Citation Index (SciSearch), Scopus, Web of Science, Chemical Abstracts, Embase, Index Copernicus, EBSCO, African Index Medicus, JournalSeek, Journal Citation Reports/Science Edition, Directory of Open Access Journals (DOAJ), African Journal Online, Bioline International, Open-J-Gate and Pharmacy Abstracts

INTRODUCTION

Orally disintegrating tablet (ODT) is an emerging and important dosage form in pharmaceutical market. It was developed by R.P. Scherer Corporation in 1986 and was first introduced into Swedish market as Zydis technology to formulate famotidine ODT in 1993[1]. Approval by FDA was done in 1996 for its use to formulate Claritin RediTabs by Schering-Plough. Currently, numerous other ODT technologies by several pharmaceutical companies and research groups are available as well [2].

Orally disintegrating tablet ODT formulations were mainly developed for existing drugs. This is because it extends product self-life, expands solid dosage form market, and avoids counterfeiting. Oral disintegrating offers ease of administration, and convenience for special

© 2023 The authors. This work is licensed under the Creative Commons Attribution 4.0 International License

populations like patients who cannot swallow tablets (pediatrics and geriatrics), patients with swallowing difficulties (dysphagia), or those who have limited access to water. Some ODT formulations allow for a high drug load, leave no grittiness or sandy feeling in mouth, allow for masking taste of bitter drugs via drug encapsulation, coating, or by using various excipients, and provide good stability. In addition convenience of administration. ODT to formulations allowed for extending product lifecycle for manufacturers, expanded market size solid dosage forms. and reduced for counterfeiting potential.

Formulation of ODT has been studied as a way of improving bioavailability of poorly watersoluble drugs and it has been observed that ODT increases the bioavailability of such drugs [3]. There are various challenges involved in formulation and development of ODT which include achieving adequate tablet hardness, accommodating high drug load, efficiently masking a bitter taste, leaving a good mouth feel following tablet disintegration, avoiding additional costs for special packaging required for friable tablets, and maintaining physical and chemical drug stability during storage [3].

Disintegration mechanisms of ODTs

Other names for ODT include orodisperse, guickdissolving system, mouth dissolving tablets, fast melting tablets, fast disintegrating tablet (FDT), fast dissolving tablet (FDT), fast dissolving tablet (RDT), rapid dissolving tablet (RDT), and orally dissolving tablet (ODT). United States Food and Drug Administration (US FDA) defined orally disintegrating tablets (ODTs) as "a solid dosage form containing medicinal substances which disintegrates rapidly, usually within a matter of seconds, when placed upon the tongue" [4]. In special patient populations like geriatrics and pediatrics, its ease of use is crucial for improving patient compliance, making it dosage form of choice [5]. Orally disintegrating tablets are made with specialized functional excipients that are meant to dissolve quickly in the mouth without need for chewing, crushing, or presence of a lot of liquid [2].

There are two main forms of ODTs namely; freeze-drving manufacturing tablets (lyophilization) and loosely compressed tablets. The freeze-drying procedure is used to increase dissolution rate and oral bioavailability of medicines that have low solubility but high permeability (biopharmaceutic classification system Class II drugs). Freeze drying (Lyophilization) involves removal of water from a product after it has been frozen. This process is carried out in a variety of ways to produce the same end product. In this method, the drug is physically trapped in a water-soluble matrix (a water-soluble mixture of saccharide and polymer formulated to provide rapid dispersion and physical strength), which is then freeze-dried to produce a product that dissolves quickly when placed in the mouth. One of the most successful freeze-drying technologies is called Zydis process. Zydis technology uses a proprietary freeze-drying technique to create completed dose units that differ greatly from standard oral systems. With this technology, a medication solution or suspension in water is poured into pre-fabricated blisters (giving tablet shape) and then frozen in a specifically engineered cryogenic freezing process to manage size of ice crystals, ensuring that the tablet has a porous matrix for quick disintegration. Frozen units are then transferred to large-scale freeze dryers for sublimation process, in which moisture from the tablets is evaporated and open blisters are sealed using heat seal method [6]. Once a tablet comes in contact with a liquid medium, the super disintegrants swell, deform, or wick up the liquid, thus increasing tablet volume. This results in break-up and disintegration. Extent of water uptake as well as rate at which water uptake develops are very important for tablets to disintegrate. Conventional equipment and limited number of steps in his method, make it most cost-effective method for manufacturing ODT. According to guidelines developed by FDA. regulation, and accreditation of ODT is an ongoing process that needs continuous review to ensure quality and efficiency for human use. Accreditation process for new drugs is usually an annual exercise though review of quality is continuous [6].

FDA ODT guidelines vs USP disintegration test

Disintegration tests are a commonly used USP method to investigate disintegration time in vivo and to assist during solid dosage form development [6]. Disintegration time is defined as time it takes for an oral dosage form to disintegrate into small particles under designed experimental conditions [7]. United States Pharmacopoeia (USP) tests are intended to provide specific quality control tests to ensure standard product performance for finished products. Therefore, current USP disintegration test using basket-rack assembly provides only pass/fail results following a specific time. In this test, a large volume of disintegration medium is utilized to evaluate tablet disintegration as well as forceful continuous vertical motion of tested tablets. This test was initially designed to test the disintegration of conventional tablets that disintegrate and release their active ingredient in the presence of large acidic aqueous volume and vigorous motility of the stomach over an extended period.

Orally disintegrating tablets (ODTs) have a disintegration time of approximately 30 s or less and recommended tablet weight should be less than 500 mg. However, USP disintegration test using basket-rack assembly does not resemble a static environment and small salivary volumes with almost a neutral pH in oral cavity. Also, due to fast disintegration of ODT formulations (< 30 s), actual disintegration end-point is easily missed by time required for ODTs to be placed in apparatus, basket to be lowered, time for apparatus to be started and stopped, and basket to be raised again for visual inspection, which definitely cannot be accounted for within 30 s specified by FDA. Furthermore, this basket-rack assembly test cannot be utilized for product development and comparison of various ODT formulations using various functional excipients and manufacturing methods [9]. Therefore, according to 2016 edition of USP, disintegration testing for ODTs is referred to as the testing utilized method and standards by the prescribed manufacturer as in individual monographs of ODT products. These guidelines have created several testing methods and standards rather than one universal method and standards for all ODTs.

Disintegration tests for ODTs

To accurately measure the disintegration time of ODTs in conditions that resemble in vivo environments, various techniques have been developed. Development of these tests was crucial for product development and optimization that facilitated studying various factors that affect rapid tablet disintegration including type and excipients used, grade of manufacturing methods, and settings. Also, it encouraged studies investigating effect of various ODT formulations containing various excipients and drua loads on pharmacokinetics and bioavailability of active ingredients. Based on factors that researchers believed to have a great impact on ODT's ability to disintegrate, different disintegration equipment designs and testing procedures were developed.

Basket-rack assembly (USP)

Basket-rack assembly is USP method used to describe ODT disintegration process and apparatus in addition to measuring an estimation

of the time it takes for an ODT to disintegrate [9,10]. Through this process, open-ended tubes are vertically held by plates and woven by stainless steel wires. Tablets are placed in tubes and lowered into beakers after which baskets are raised and lowered into 900 mL immersion medium at constant frequencies until full disintegration. The process is repeated on tablets that take long to disintegrate to help calculate the average time it takes for such tablets to fully disintegrate. This method is good as it allows for physical observation of disintegration process in addition to allowing multiple experimental performances help determine to actual disintegration time [9].

Basket-rack apparatus is prepared in a uniformly standard measure to help provide consistent results. This apparatus may be varied in design but specifications of glass tubes and screen meshes need to be constantly maintained [9]. Disintegration process is monitored through a glass apparatus and the end-time of operation is recorded. This method has the challenge of being a confusing process, has no good correlation with *in vivo* tablet disintegration, and volume of medium used is not applicable in oral cavity.

Novel method utilizing a CCD camera

This new method involves use of CCD cameras to help evaluate disintegration of ODTs. It was designed to help complement the Japanese Pharmacopoeia (JP) XIII disintegration test that measures the dissolution of orally disintegrating tablets among elderly people in society [11]. This study is predicated on observation that methods for detecting disintegration rate of ODTs, such as JP XII, are ineffective due to agitation impact they produce. As a result, this method complementing JP XIII was developed to help ensure an efficient measure of disintegration time [11]. In this method, a camera is interfaced along a path of disintegration with a personal computer and a motion sensor to help monitor all movements. Thereafter, CCD camera records disintegration course while transferring pictures to a computer for analysis and calculation of actual disintegration time [11].

The novel CCD camera is a technique considered an advancement from other *in vivo* and *in vitro* measurement methods. It takes into account various components of ODTs like binders, and parameters such as thickness and diameter to determine exact disintegration time [11]. Disintegration process takes place in a slightly agitated medium giving room for disintegration of RDTs of different formulations.

This technique is very advanced and efficient as it takes into account morphological changes occurring during the disintegration process [12]. This suggests why it is appropriately used in quality control of ODT production process and other pharmaceutical development processes [11].

In another innovative disintegration test method, Kenada et al. [13] also included use of a CCD camera to monitor disintegration time of RDTs. In this method, CCD camera was used to monitor behavioral changes in disintegrating tablets. This helps monitor disintegrating test also mechanisms in addition to accurately determining exact disintegration time as it precisely records start and end of reaction time. Despite these benefits, it is difficult to construct an efficient CCD camera accommodating apparatus.

Texture analyzer (TA)

Through a TA process, a new operating structure is developed and designed like a human oral system through which attached granules are systematically eliminated during disintegration It is a process of measuring process. disintegration time of RDTs in an artificial environment similar to oral cavity [14]. In texture analyzer method, RDT is placed on a grid and lowered into a disintegration medium by a trigger force that pushes it down a disintegration medium [14]. Effect of changing medium and temperature is assessed during evaluation. An increase in temperature and nature of medium affect disintegration time of RDTs as observed by use of artificial saliva and distilled water mediums [14]. Quantitative values are calculated through disintegration profile of RDTs with measurements of starting time being recorded as T1 and end process time recorded as T₂. Values gained through TA aid in obtaining RDT profiles in addition to identifying an excellent correlation with *in-vivo* disintegration time in a patient's mouth. This process is largely applauded for its high accuracy level that enables it to provide a mouth-feel of RDT as it establishes thickness of tablets with disintegration distance obtained from their profiles [14]. It is also advantageous as it allows for differentiation of various rapidly disintegrating tablets through their various characteristics as outlined in their profiles. However, it is difficult to always construct exact accurate equipment with standard measurements. So, disintegration time of these tablets using TA is very short, making it difficult to give an accurate measure of disintegration time using a stopwatch.

The Novel disintegration apparatus (Rotary Shaft)

It is simpler to ascertain disintegration time of various RDT formulations using this disintegration-time device [15]. This approach is faster compared to other conventional approaches.

The method involves use of a rotating shaft that has a weight of about 10 g and a speed of rotation of about 10 rpm. It employs use of tablets stored in both severe conditions and appropriate conditions of about 37 °C Disintegration time of tablets obtained using this method provides values similar to those of human sensory tests conducted in the oral cavity involving saliva and tongue pressure towards the mouth palate [15]. upper This test is advantageous because it provides accurate results that directly correlate to what exactly happens to RDTs in human oral cavity [16]. The approach is difficult to use for precise observations because it is impossible to witness the mechanism that causes pill disintegration. The same test was adopted in another study in which direct compression orally disintegrating tablets were applied in an apparatus set that provides different volumes of purified water as dissolution medium in different temperatures [17]. However, the provision of magnetic resonance imaging (MRI) helps in providing an analysis of the properties of disintegrating tablets, helps to determine exact disintegration time, and improves quality of tested results as well as mechanism of disintegration process. With this apparatus, it is easier to determine effect of compression force, temperature, and volume of medium in enhancing disintegration [17]. This method has the challenge of requiring extremely high attention levels.

Prototypical disintegration tester (PDT)

The prototypical disintegration tester (PDT) measures in-vitro disintegration time for ODTs with high accuracy [13]. To help achieve more advanced results, this disintegration method employs use of a CCD camera to monitor behaviour and reaction mechanisms of ODT profiles. Furthermore, CCD camera attached to this method helps enhance accurate determination of reaction end-time. With this process, water absorption rate of ODTs is measured in an artificial environmental set-up with a commercial tester, and an estimation of disintegration time was measured. Using swelling and water absorption rate, this method is used to accurately determine disintegration time of ODT with an additional enhancement of a CCD

camera [13]. The method is particularly helpful since it provides a precise estimation of disintegration time and aids in describing disintegrating behaviour. Difficulty in creating the test apparatus is a major draw-back [13].

Novel method using commercial RDTs

The novel method of measuring disintegration time is an advancement of conventional method described in Japanese Pharmacopoeia (JP) XIII [18]. Many studies were done to determine disintegration time associated with rapidly disintegrating tablets, in vitro using this new and while comparing it with simple method conventional disintegration test in JP XII which is inefficient due to its environmental conditions that largely differ from those in oral cavity, and with human sensory test [19]. New method involves use of suspension meshes kept at room temperature sandwiching the RDT with some additional weight exerted to facilitate disintegration process. A drop of test fluid, which is either purified water or artificial saliva induces wetting for disintegration to commence. This process equally takes a relatively short time, comparable to that taken by healthy volunteers, which means that this new method is most likely efficient and reliable [18].

Conventional method on the other hand utilizes a large amount of water that is not present in oral cavity with no associated pressure associated resulting in a slow disintegration time. Novel method incorporates main factors that facilitate disintegration, which include volume of oral saliva and internal pressure that pushes against the mouth palate [18]. Thus, novel method offers special benefits such as ease of construction, adjustable oral cavity environment, and ease of measuring at least two aspects of tablets. Furthermore, medium in this method impinges on the tablet drop by drop and also only from topside, which makes it unreliable and does not reflect *in vivo* disintegration time.

OD-mate testing

The OD-mate has been employed to obtain disintegration time during extraction of amlodipine ODT formulations [20]. Orallv disintegrating tablet (ODT) was placed on a trapezoidal mesh in a flat-bottomed test tube and compressed by two weights [21]. Measurement commences as soon as water enters the test tube, and time to completely disintegrate is recorded [20]. In a study, six (6) women who had not experienced any respiratory issues in week before test were subjected to a gustatory test to

determine how bitter different quinine hydrochloride concentrations tasted in their predetermined seconds mouths at and disintegration obtained [20]. Before test, subjects were instructed to thoroughly rinse their mouths [22]. Disintegration time was measured using a chronometer since it is felt when ODT granules are completely dissolved in buccal cavity [20]. Test subjects were given time to gargle thoroughly and wait for about 20 mins before they could take another test, which determined how accurate the approach was. Pressing of granules against upper palate surface of the mouth was encouraged, though swallowing of period saliva durina test was stronalv discouraged. Disintegration rate was evaluated by using both OD-mate and gustatory sensation while the taste of product was believed to estimate the bitterness intensity of amlodipine ODT formulations [23,24]. The method is very accurate as evidenced by strong correlation of results between artificial taste sensors and gustatory sensation testing techniques. However, the weight placed on tablets (which may be too high) and the absence of a rotary shaft or pressure to enhance disintegration process are major drawbacks [20].

Modified wetting test

Wetting test is a commonly used method of assessing disintegration time of orally disintegrating tablets (ODTs) [8]. This test involves immersing samples of ODTs in an appropriate amount of water and recording disintegration time. Amount of water used in wetting test experiment is assumed equivalent to that in human oral cavity. In a previous analysis, most studies noted that wetting test experiments failed to show exact oral cavity environment [8]. Water absorption ratio while determining disintegration rate was calculated to improve accuracy of these results. However, another study found that even though certain tablets will swell after absorbing water and fail to disintegrate, it may be possible to estimate disintegration time by wetting test [8]. This, therefore, was a major challenge of wetting test experiment, hence, allowing for use of a more effective and innovative method, which is the shaking bath test method [8]. With this method, all types of ODT tablets are effectively measured and exact disintegration time are recorded. Enhanced shaking test, which improved upon the previous wetting test method, allowed for measurement of disintegration times rather than water permeation times. The only difficulty is that water absorption ratios are not measured.

Bafail

Table 1: Methods for measuring disintegration time

Method	Advantages	Disadvantages
Basket rack assembly (USP)	It allows for physical observation of the disintegration process in addition to allowing multiple experimental performances to help determine actual disintegration pace.	It is a confused process.
Novel method utilizi ng CCD camera	Outlines RDTs morphological changes during disintegration process and enhances accuracy and quality of monitoring disintegration time.	It is difficult to construct an efficient CCD camera accommodating apparatus.
Texture analyzer (TA)	It allows for differentiation of various RDTs.	Complexity to construct an apparatus with exact standard measurements.
Novel disintegration apparatus	It is accurate in measuring disintegration time, ability to determine effect of compression force, temperature, and volume of medium enhancing disintegration process	Requires a lot of attention.
Prototypical disintegration tester (PDT)	It provides accurate results that are directly correlated to <i>in vivo</i> disintegration in oral cavity and, also ability to determine disintegrating mechanisms.	Developing test apparatus is considered a complex process.
Novel and simple method	Easier to construct, adjustable oral cavity environment and also measures various aspects of ODTs	Disintegration time strongly depends on available oral cavity environment and conditions that are difficult to measure.
OD-Mate Testing	It is very accurate as evidenced by strong correlation of results between artificial taste sensors and gustatory sensation testing techniques.	Weight above tablet is very high and there is no rotary shaft or pressure to enhance disintegration process.
Modified Wetting Test	Disintegration time of all ODTs is effectively measured	Inability to measure water absorption ratio of ODTs.

Standardized wetness shake tests are highly effective in determining times at which different ODTs disintegrate, providing more accurate results [8].

The methods presented with their characteristic advantages that enable efficient estimation of disintegration time are shown in Table 1.

CONCLUDING REMARKS

Introduction of ODTs has solved some problems associated with drug administration to pediatric and elderly patients, who account for a sizable proportion of world's population. Number of companies in ODT drug delivery market is high, as evidenced by an ever increasing number of patents and marketed products. Among other drug delivery companies, those in ODT sector have a huge potential for extending product selflife and increasing profitability of existing pharmaceuticals. This delivery facilitates a very short disintegration time. Because disintegration time is an essential parameter for ODT development, regular USP disintegration test is not suitable. All of the studies developed are evidence that there is a variety of effective methods to evaluate disintegration time and mechanisms. It is therefore an expectation that future trends in drug delivery innovation will continue to bring together different disciplines to create novel technologies.

DECLARATIONS

Acknowledgements

None provided.

Funding

None provided.

Ethical approval

None provided.

Availability of data and materials

The datasets used and/or analyzed during the current study are available from the corresponding author on reasonable request.

Conflict of Interest

No conflict of interest associated with this work.

Contribution of Authors

The authors declare that this work was done by the authors named in this article and all liabilities pertaining to claims relating to the content of this article will be borne by them.

Open Access

This is an Open Access article that uses a funding model which does not charge readers or their institutions for access and distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/ 4.0) and the Budapest Open Access Initiative (http://www.budapestopenaccessinitiative.org/rea d), which permit unrestricted use, distribution, and reproduction in any medium, provided the original work is properly credited.

REFERENCES

- Rathbone MJ, Hadgraft J, Roberts MS. "The Zydis Oral Fast-Dissolving Dosage Form" Modified-Release Drug Delivery Technology. Modified-Release Drug Delivery Technology CRC Press. 2002; p 191
- 2. Vinushitha VSaS. Oral Disintegrating Tablets: An Overview. Int J Chem Pharm Sci 2010; 1(2): 1-12.
- Nayak AK, Manna K. Current developments in orally disintegrating tablet technology. J Pharm Educ Res 2011; 2(1): 21-34
- U.S. Department of Health and Human Services: Food and Drug Administration. Center for Drug Evaluation and Research (CDER). (2008). "Guidance for Industry Orally Disintegrating Tablets." P. 2.
- 5. Klancke J. Dissolution Testing of Orally Disintegrating Tablets. Dissolution Technol 2003; 10(2): 6-9.
- Shoukri RA, Ahmed IS, Shamma RN. In vitro and in vivo evaluation of nimesulide lyophilized orally disintegrating tablets. Eur J Pharm Biopharm: Official Journal of Arbeitsgemeinschaft fur Pharmazeutische Verfahrenstechnik eV 2009; 73(1): 162-171.
- Prajapati GB, Patel B. Formulation, Evaluation and Optimization of Orally Disintegrating Tablet of Piroxicam. Int J Pharm Tech Res 2010; 2(2): 1893-1899.
- 8. Hooper P, Lasher J, Alexander KS, Baki G. A new modified wetting test and an alternative disintegration

test for orally disintegrating tablets. J Pharm Biomed Anal 2016; 120: 391-396.

- USB/NF. General Chapters: Disintegration. United States Pharmacopeia. 1. 37/32 ed. William E. Brown: United States Pharmacopeial Convention Ip.
- Que L, Wu W, Cheng X, Hu T. Evaluation of Disintegrating Time of Rapidly Disintegrating Tablets by a Paddle Method. Pharm Dev Technol 2006; 11(3): 295-301.
- Morita Y, Tsushima Y, Yasui M, Termoz R, Ajioka J, Takayama K. Evaluation of the disintegration time of rapidly disintegrating tablets via a novel method utilizing a CCD camera. Chem Pharm Bull 2002; 50(9): 1181-1186.
- Khan S, Kataria P, Nakhat P, Yeole P. Taste masking of ondansetron hydrochloride by polymer carrier system and formulation of rapid-disintegrating tablets. AAPS PharmSciTech 2007; 8(2): Article 46.
- Kondo K, Niwa T, Danjo K. Evaluation of disintegration properties of orally rapidly disintegrating tablets using a novel disintegration tester. Chem Pharm Bull 2012; 60(10): 1240-1248.
- Abdelbary G, Eouani C, Prinderre P, Joachim J, Reynier J, Piccerelle P. Determination of the in vitro disintegration profile of rapidly disintegrating tablets and correlation with oral disintegration. Int J Pharma 2005; 292(1-2): 29-41.
- 15. Harada T, Narazaki R, Nagira S, Ohwaki T, Aoki S, Iwamoto K. Evaluation of the disintegration properties of commercial famotidine 20 mg orally disintegrating tablets using a simple new test and human sensory test. Chem Pharm Bull 2006; 54(8): 1072-1075.
- 16. Garbacz G, Wedemeyer RS, Nagel S, Giessmann T, Monnikes H, Wilson CG, et al. Irregular absorption profiles observed from diclofenac extended-release tablets can be predicted using a dissolution test apparatus that mimics in vivo physical stresses. Eur J Pharma Biopharma: Official journal of Arbeitsgemeinschaft fur Pharmazeutische Verfahrenstechnik eV 2008; 70(2): 421-428.
- Brniak W, Jachowicz R, Krupa A, Skorka T, Niwinski K. Evaluation of co-processed excipients used for direct compression of orally disintegrating tablets (ODT) using novel disintegration apparatus. Pharm Dev Technol 2013; 18(2): 464-474.
- Hoashi Y, Tozuka Y, Takeuchi H. Development of a novel and simple method to evaluate disintegration of rapidly disintegrating tablets. Chem Pharma Bull 2013; 61(9): 962-926.
- Kakutani R, Muro H, Makino T. Development of a new disintegration method for orally disintegrating tablets. Chem Pharm Bull 2010; 58(7): 885-890.
- 20. Uchida T, Yoshida M, Hazekawa M, Haraguchi T, Furuno H, Teraoka M, Ikezaki H. Evaluation of palatability of 10 commercial amlodipine orally disintegrating tablets by gustatory sensation testing, OD-mate as a new disintegration apparatus, and the artificial taste sensor. J Pharm Pharmacol 2013; 65(9): 1312-1320.

Trop J Pharm Res, November 2023; 22(11): 2405

- Yoshida M, Haraguchi T, Uchida T. Bitterness evaluation of acidic pharmaceutical substances (NSAIDs) using a taste sensor. Chem Pharm Bull 2014; 62(12): 1252-1258.
- 22. Yano K, Mita S, Morimoto K, Haraguchi T, Arakawa H, Yoshida M, Yamashita F, Uchida T, Ogihara T. Multiple Linear Regression Analysis Indicates Association of P-Glycoprotein Substrate or Inhibitor Character with Bitterness Intensity, Measured with a Sensor. J Pharm Sci 2015; 104(9): 2789-2794.
- Anand V, Kataria M, Kukkar V, Saharan V, Choudhury PK. The latest trends in the taste assessment of pharmaceuticals. Drug Discov Today 2007; 12(5–6): 257-265.
- Tokuyama E, Matsunaga C, Yoshida K, Mifsud JC, Irie T, Yoshida M, Uchida T. Famotidine orally disintegrating tablets: bitterness comparison of original and generic products. Chem Pharm Bull 2009; 57(4): 382-387.
- Rawas-Qalaji MM, Simons FE, Simons KJ. Fastdisintegrating sublingual tablets: effect of epinephrine load on tablet characteristics. AAPS PharmSciTech 2006; 7(2): E41.
- Narang N, Sharma J. Sublingual mucosa as a route for systemic drug delivery. Int J Pharm Pharm Sci 2011; 3(2): 18-22.

- Narazaki R, Harada T, Takami N, Kato Y, Ohwaki T. A new method for disintegration studies of rapid disintegrating tablet. Chem Pharm Bull 2004; 52(6): 704-707.
- Kraemer J, Gajendran J, Guillot A, Schichtel J, Tuereli A. Dissolution testing of orally disintegrating tablets. J Pharm Pharmacol 2012; 64(7): 911-918.
- Ozeki Y, Watanabe Y, Inoue S, Danjo K. Comparison of the compression characteristics between new one-step dry-coated tablets (OSDRC) and dry-coated tablets (DC). Int J Pharm 2003; 259(1-2): 69-77.
- 30. Miyazaki S, Kubo W, Itoh K, Konno Y, Fujiwara M, Dairaku M, Togashi M, Mikami R, Attwood D. The effect of taste masking agents on in situ gelling pectin formulations for oral sustained delivery of paracetamol and ambroxol. Int J Pharm 2005; 297(1-2): 38-49.
- Haraguchi T, Yoshida M, Uchida T. Evaluation of ebastine-loaded orally disintegrating tablets using new apparatus of detecting disintegration time and e-tongue system. J Drug Deliv Sci Technol 2014; 24(6): 684-688.
- 32. Bi Y, Sunada H, Yonezawa Y, Danjo K, Otsuka A, Iida K. Preparation and evaluation of a compressed tablet rapidly disintegrating in the oral cavity. Chem Pharm Bull 1996; 44(11): 2121-2127.